

ZYKLUSZEIT

00:34:28

H : MIN : SEK

STERILISATIONSDAUER

3:30

MIN : SEK

STERILISATION TEMP.

134.1 134.8

MIN MAX C

F0

12.5

LETZTE TESTS

BOWIE-DICK LEAK TEST

RACK-ID

-

PROGRAMM-NR. / VERSION

-

PROGRAMM-NAME

Bowie-Dick-Test

PROZESSDATEN (Zusammenfassung)

Zeit	Phase	T1 C	T2 C	P1 mbar	P2 mbar
0:02	1. Vorvakuum	46.0	46.0	1034	
3:00	1. Druck	55.9	55.9	2050	
5:00	2. Vorvakuum	65.2	65.2	55	
7:30	2. Druck	80.1	80.1	2050	
10:00	3. Vorvakuum	95.4	95.4	70	
12:30	3. Druck	105.8	105.8	2050	
15:00	4. Vorvakuum	110.3	110.3	65	
17:30	4. Druck	120.5	120.5	2050	
20:00	Aufheizung	130.2	130.2	1850	
22:00	Sterilisation	134.1	134.1	2100	
31:15	Druckausgleich	62.3	62.3	984	
33:00	Trocknung	68.1	68.1	85	
34:30	Programm Ende	69.2	69.2	984	

OPERATOR

User

CHARGE

freigegeben

☐

ja

☐

nein

☐

bedingte Freigabe

BEMERKUNGEN ZUR BEDINGTEN FREIGABE

SIGNATURE

Zeit	Meld. Nr.	Meldung	Phase	T1 C	T2 C	P1 bar a	P2 bar a	P3 bar a	P4 bar a
0:02			1. Vorvakuum	46.00	46.00	1034	1034		
3:00			1. Druck	55.90	55.90	2050	2050		
5:00			2. Vorvakuum	65.20	65.20	55	55		
7:30			2. Druck	80.10	80.10	2050	2050		
10:00			3. Vorvakuum	95.40	95.40	70	70		
12:30			3. Druck	105.80	105.80	2050	2050		
15:00			4. Vorvakuum	110.30	110.30	65	65		
17:30			4. Druck	120.50	120.50	2050	2050		
20:00			Aufheizung	130.20	130.20	1850	1850		
22:00			Sterilisation	134.10	134.10	2100	2100		
25:00			Sterilisation	134.50	134.50	2150	2150		
28:30			Sterilisation	134.80	134.80	2180	2180		
31:15			Druckausgleich	62.30	62.30	984	984		
33:00			Trocknung	68.10	68.10	85	85		
34:30			Programm Ende	69.20	69.20	984	984		